

Anti-KRT86 Antibody

Catalog Number: A08530-1

About KRT86

KRT86 is a type II keratin protein, which heterodimerizes with type I keratins to form hair and nails. This gene is present in a cluster of related genes and pseudogenes on chromosome 12. Mutations in this gene have been observed in patients with the hair disease monilethrix.

Overview

Product Name	Anti-KRT86 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-KRT86 Antibody catalog # A08530-1. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O43790

Technical Details

Immunogen	E.coli-derived human KRT86 recombinant protein (Position: A442-C483).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-KRT86 Antibody

For Research Use Only. Not for use in diagnostic procedures.