

Anti-BRS-3 (C203) Antibody

Catalog Number: A08571-2


Overview

Product Name	Anti-BRS-3 (C203) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-BRS-3 (C203) Antibody catalog # A08571-2. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P32247

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 280-330 of Human BRS-3.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IF: 1:50-1:200

Anti-BRS-3 (C203) Antibody (A08571-2) Images



Western blot (WB) analysis of BRS-3 (C203) pAb at 1:500 dilution Lane1:A549 whole cell lysate(10ug) Lane2:Panc1 whole cell lysate(10ug) Lane3:BV2 whole cell lysate(40ug) Lane4:H9C2 whole cell lysate(40ug) Lane5:PC12 whole cell lysate(40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-BRS-3 (C203) Antibody

For Research Use Only. Not for use in diagnostic procedures.