

Anti-GAS2 Antibody

Catalog Number: A08589

About GAS2

Plays a critical role in endothelial cell capillary morphogenesis By similarity.

St Croix B., Science 289:1197-1202(2000).

Carson-Walter E.B., Cancer Res. 61:6649-6655(2001).

Nanda A., Cancer Res. 64:8507-8511(2004)

Overview

Product Name	Anti-GAS2 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-GAS2 Antibody catalog # A08589. Tested in WB,IHC applications. This antibody reacts with Human,Mouse.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O43903

Technical Details

Immunogen	Synthesized peptide derived from Internal of human PLXDC1.
Predicted Reactive Species	Chimpanzee, Drosophila, Macaque
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.

If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.
Some PubMed article(s) citing the expression level of this target are as follows:
Boster Bio's internal QC testing used:
WB: 1:500-1:2000
IHC: 1:50-1:200

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GAS2 Antibody