

Anti-RPL10A Antibody

Catalog Number: A08606-1

Overview

| | |
|----------------------|--|
| Product Name | Anti-RPL10A Antibody |
| Reactive Species | Human, Mouse |
| Description | Boster Bio Anti-RPL10A Antibody catalog # A08606-1. Tested in WB, ICC, IF applications. This antibody reacts with Human, Mouse. |
| Application | IF, ICC, WB |
| Clonality | Polyclonal |
| Formulation | 500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | 12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing. |
| Host | Rabbit |
| Uniprot ID | P62906 |

Technical Details

| | |
|---------------------|--|
| Immunogen | E.coli-derived human RPL10A recombinant protein (Position: S3-Q214). |
| Form | Liquid |
| Concentration | 500 ug/ml |
| Purification | Immunogen affinity purified. |
| Suggested Dilutions | Western blot, 1:500-2000 Immunocytochemistry/Immunofluorescence, 1:50-400 |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RPL10A Antibody

For Research Use Only. Not for use in diagnostic procedures.