

Anti-HMG20B Antibody

Catalog Number: A08646-3

About HMG20B

CoREST (corepressor of repressor element-1-silencing transcription factor) is a chromatin- modifying complex that was first identified as a regulator of neuronal gene expression, having a essential role in neuronal differentiation [PMID:10449787]. The CoREST complex is recruited by REST (RE1-silencing transcription factor) that recognizes the 21-nucleotide RE-1 sequence in the promoter of more than 1,700 different target genes in human and mouse [PMID:15240883]. HMG20B is one component of the CoREST complex, which mediates chromatin remodeling and repression of specific genes involved in neuronal differentiation through modulation of chromatin structure [PMID:19997604].

Overview

Product Name	Anti-HMG20B Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-HMG20B Antibody catalog # A08646-3. Tested in WB, IHC, ICC, IF, IP, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9P0W2

Technical Details

Immunogen	E.coli-derived human HMG20B recombinant protein (Position: T75-F282).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400

ImmunoPrecipitation, 1:250-300
ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HMG20B Antibody

For Research Use Only. Not for use in diagnostic procedures.