

Anti-Leiomodin-1 LMOD1 Antibody

Catalog Number: A08666


Overview

Product Name	Anti-Leiomodin-1 LMOD1 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Leiomodin-1 LMOD1 Antibody catalog # A08666. Tested in WB,IHC applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P29536

Technical Details

Immunogen	Recombinant protein of human LMOD1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC 1:50-1:200

Anti-Leiomodin-1 LMOD1 Antibody (A08666) Images



Immunohistochemistry of paraffin-embedded human colon using LMOD1 antibody at dilution of 1:200 (40x lens).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Leiomodin-1 LMOD1 Antibody

For Research Use Only. Not for use in diagnostic procedures.