

Anti-Arylsulfatase D ARSD Antibody

Catalog Number: A08673-1

Overview

Product Name	Anti-Arylsulfatase D ARSD Antibody
Reactive Species	Human
Description	Boster Bio Anti-Arylsulfatase D ARSD Antibody catalog # A08673-1. Tested in ELISA, IF, IHC applications. This antibody reacts with Human.
Application	ELISA, IF, IHC
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P51689

Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human ARSD. AA range:331-380
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	IHC 1:100-1:300 IF 1:200-1:1000 ELISA 1:20000

Anti-Arylsulfatase D ARSD Antibody (A08673-1) Images




Figure 1. Immunohistochemistry validation of ARSD using Anti-Arylsulfatase D ARSD Antibody (A08673-1).

Immunohistochemical analysis of paraffin-embedded human lung cancer. Antibody was diluted at 1:100 (4°C)

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Arylsulfatase D ARSD Antibody