

Anti-Calponin 2 (D156) CNN2 Antibody

Catalog Number: A08680-1


Overview

Product Name	Anti-Calponin 2 (D156) CNN2 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Calponin 2 (D156) CNN2 Antibody catalog # A08680-1. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q99439


Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 122-170 of Human Calponin 2.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IF: 1:50-1:200

Anti-Calponin 2 (D156) CNN2 Antibody (A08680-1) Images



Western blot (WB) analysis of Calponin 2 (D156) pAb at 1:500 dilution Lane1:AML-12 whole cell lysate(40ug) Lane2:H9C2 whole cell lysate(20ug) Lane3:L02 whole cell lysate(40ug) Lane4:HepG2 whole cell lysate(40ug) Lane5:Hela whole cell lysate(40ug)



Immunohistochemistry (IHC) analyzes of Calponin 2 (D156) pAb in paraffin-embedded human brain tissue.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Calponin 2 (D156) CNN2 Antibody

For Research Use Only. Not for use in diagnostic procedures.