

## Anti-DRAK2 (I286) STK17B Antibody

Catalog Number: A08711-3

### Overview


|                      |  |
|----------------------|--|
| Product Name         | Anti-DRAK2 (I286) STK17B Antibody  |
| Reactive Species     | Human, Mouse, Rat  |
| Description          | Boster Bio Anti-DRAK2 (I286) STK17B Antibody catalog # A08711-3. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat. |
| Application          | IHC, WB  |
| Clonality            | Polyclonal   |
| Formulation          | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2   |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host                 | Rabbit   |
| Uniprot ID           | O94768   |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | Synthetic peptide, corresponding to amino acids 261-310 of Human DRAK2.   |
| Isotype             | IgG   |
| Form                | Liquid  |
| Concentration       | 1 mg/ml   |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | WB: 1:500-1:1000<br>IHC: 1:50-1:200   |

## Anti-DRAK2 (I286) STK17B Antibody (A08711-3) Images

---



Immunohistochemistry (IHC) analyzes of DRAK2 (I286) pAb in paraffin-embedded human breast carcinoma tissue at 1:50.

### Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DRAK2 (I286) STK17B Antibody

For Research Use Only. Not for use in diagnostic procedures.