

Anti-Coatomer subunit delta ARCN1 Antibody

Catalog Number: A08725

About ARCN1

Phosphorylates myosin light chains By similarity. Acts as a positive regulator of apoptosis.

Sanjo H., J. Biol. Chem. 273:29066-29071(1998).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Oppermann F.S., Mol. Cell. Proteomics 8:1751-1764(2009).

Overview

Product Name	Anti-Coatomer subunit delta ARCN1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Coatomer subunit delta ARCN1 Antibody catalog # A08725. Tested in WB,ICC/IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P48444

Technical Details

Immunogen	Synthesized peptide derived from internal of human STK17B.
Predicted Reactive Species	Chimpanzee, Drosophila, Macaque
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.

If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.

Some PubMed article(s) citing the expression level of this target are as follows:

Boster Bio's internal QC testing used:

WB: 1:500-1:2000

ICC/IF: 1:50-1:200

Anti-Coatomer subunit delta ARCN1 Antibody (A08725) Images




Figure 1. Western blotting validation for Anti-Coatomer subunit delta ARCN1 Antibody A08725

Western blot analysis of extracts of various cell lines

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Coatomer subunit delta ARCN1 Antibody