

## Anti-TS1R2 TAS1R2 Antibody

Catalog Number: A08729

### Overview


Product Name	Anti-TS1R2 TAS1R2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-TS1R2 TAS1R2 Antibody catalog # A08729. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8TE23

### Technical Details

Immunogen	Synthesized peptide derived from human protein at AA range: 190-270
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	TS1R2 antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-2000 ELISA 1:5000-20000

## Anti-TS1R2 TAS1R2 Antibody (A08729) Images

---



Western Blot analysis of Hela lysis, using primary antibody at 1:1000 dilution. Secondary antibody was diluted at 1:10000

---

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TS1R2 TAS1R2 Antibody

For Research Use Only. Not for use in diagnostic procedures.