

Anti-Cystatin-SA CST2 Antibody

Catalog Number: A08736-1


Overview

| | |
|----------------------|--|
| Product Name | Anti-Cystatin-SA CST2 Antibody |
| Reactive Species | Human, Mouse |
| Description | Boster Bio Anti-Cystatin-SA CST2 Antibody catalog # A08736-1. Tested in WB applications. This antibody reacts with Human,Mouse. |
| Application | WB |
| Clonality | Polyclonal |
| Formulation | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | P09228 |

Technical Details

| | |
|---------------------|---|
| Immunogen | Recombinant fusion protein of human CST2(NP_001313.1). |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | WB: 1:500-1:2000 |

Anti-Cystatin-SA CST2 Antibody (A08736-1) Images



Western blot analysis of extracts of mouse heart, using CST2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 5min.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cystatin-SA CST2 Antibody

For Research Use Only. Not for use in diagnostic procedures.