

## Anti-FoxL1 (G150) Antibody

Catalog Number: A08745

### Overview


Product Name	Anti-FoxL1 (G150) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-FoxL1 (G150) Antibody catalog # A08745. Tested in IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q12952

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 120-170 of Human FoxL1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200

## Anti-FoxL1 (G150) Antibody (A08745) Images

---



Immunohistochemistry (IHC) analyzes of FoxL1 (G150) pAb in paraffin-embedded human lung carcinoma tissue at 1:50, showing Cytoplasm and Nuclear staining. Negative control (the right) Using PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG-biotin followed by avidin-peroxidase.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FoxL1 (G150) Antibody

For Research Use Only. Not for use in diagnostic procedures.