

## Anti-SLC22A9 Antibody Picoband® Cy3 Conjugated

Catalog Number: A08762-1-Cy3

### About SLC22A9

Solute carrier family 22 member 9 is a protein that in humans is encoded by the SLC22A9 gene. Enables anion:anion antiporter activity; short-chain fatty acid transmembrane transporter activity; and sodium-independent organic anion transmembrane transporter activity. Involved in hormone transport; short-chain fatty acid import; and sodium-independent organic anion transport. Located in basolateral plasma membrane.

### Overview

Product Name	Anti-SLC22A9 Antibody Picoband® Cy3 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q8IVM8

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human SLC22A9.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SLC22A9 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.