

## Anti-FAM83D Antibody Picoband® FITC Conjugated

Catalog Number: A08800-1-FITC

### About FAM83D

Enables kinesin binding activity; microtubule binding activity; and protein kinase binding activity. Involved in several processes, including positive regulation of cell cycle G1/S phase transition; protein localization to mitotic spindle; and regulation of intracellular signal transduction. Located in cytosol; intercellular bridge; and mitotic spindle pole.

### Overview

Product Name	Anti-FAM83D Antibody Picoband® FITC Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9H4H8

### Technical Details

Immunogen	E.coli-derived human FAM83D recombinant protein (Position: P44-H559).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-FAM83D Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.