

Anti-SEZ6L Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A08803-1-Fluoro488

About SEZ6L

Seizure 6-like protein is a protein that in humans is encoded by the SEZ6L gene. SEZ6L (Seizure 6-Like protein; also BSRP-B in rodent) is a 200 kDa member of the SEZ6 family of proteins. It is expressed in the plasma membrane and ER of CNS neurons and lung bronchiolar epithelium, and likely regulates PKC alpha activation. Mature human SEZ6L is a 996 amino acid (aa) type I transmembrane protein. It contains a 930 aa extracellular domain (ECD) (aa 29-958) plus a 45 aa cytoplasmic tail. The ECD shows three CUB domains and five short consensus repeats (SCRs). There are multiple isoforms possible for SEZ6L. One shows a deletion of aa 933-943, and a second shows the same deletion coupled to a Thr substitution for aa 803-867. A third isoform contains a Leu substitution for aa 867-943. An alternative start site is reported at Met228. This start site may occur individually, or be accompanied by a deletion of aa 932-944, the same deletion of aa 932-944 coupled to a Thr substitution for aa 803-867, or a dual deletion of aa 869-945 and 1016-1024. Over aa 29-958, human BSRP-B is 80% aa identical to mouse BSRP-B.

Overview

Product Name	Anti-SEZ6L Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9BYH1

Technical Details

Immunogen	E.coli-derived human SEZ6L recombinant protein (Position: E662-E953).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488

	Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SEZ6L Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.