

Anti-HSPB11/IFT25 Antibody

Catalog Number: A08810-2

About IFT25

Predicted to enable metal ion binding activity. Predicted to be involved in kidney development and spermatogenesis. Predicted to act upstream of or within animal organ development; left/right axis specification; and skeletal system development. Predicted to be located in ciliary tip. Predicted to be part of intraciliary transport particle B. Predicted to be active in centrosome and cilium.

Overview

Product Name	Anti-HSPB11/IFT25 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-HSPB11/IFT25 Antibody catalog # A08810-2. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9Y547

Technical Details

Immunogen	E.coli-derived human HSPB11 recombinant protein (Position: M1-S140).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HSPB11/IFT25 Antibody

For Research Use Only. Not for use in diagnostic procedures.