

## Anti-POLM Antibody Picoband® HRP Conjugated

Catalog Number: A08812-2-HRP

### About POLM

Predicted to enable DNA-directed DNA polymerase activity. Predicted to be involved in double-strand break repair via nonhomologous end joining. Predicted to act upstream of or within B cell differentiation and somatic hypermutation of immunoglobulin genes. Predicted to be located in nucleoplasm. Predicted to be active in nucleus.

### Overview

Product Name	Anti-POLM Antibody Picoband® HRP Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9NP87

### Technical Details

Immunogen	E.coli-derived human POLM recombinant protein (Position: E36-H481). Human POLM shares 76.5% amino acid (aa) sequence identity with mouse POLM.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-POLM Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.