

## Anti-EphA6/Eph Receptor A6 Monoclonal Antibody

Catalog Number: A08818


### Overview

Product Name	Anti-EphA6/Eph Receptor A6 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-EphA6/Eph Receptor A6 Monoclonal Antibody catalog # A08818. Tested in WB, ELISA applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Monoclonal 3D5B1
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q9UF33

### Technical Details

Immunogen	Purified recombinant fragment of EphA6 (aa695-795) expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Ascitic fluid
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

## Anti-EphA6/Eph Receptor A6 Monoclonal Antibody (A08818) Images



Western Blot (WB) analysis using EphA6 Monoclonal antibody against truncated MBP-EphA6 recombinant protein (1) and truncated GST-EphA6(aa695-795) recombinant protein (2).

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-EphA6/Eph Receptor A6 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.