

Anti-PLB1 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A08822-1-Fluoro488

About PLB1

Predicted to enable carboxylic ester hydrolase activity. Acts upstream of or within positive regulation of acrosome reaction. Predicted to be located in apical plasma membrane. Predicted to be active in brush border membrane. Is expressed in diencephalon lateral wall mantle layer; hypothalamus mantle layer; lung; midbrain mantle layer; and telencephalon mantle layer. Orthologous to human PLB1 (phospholipase B1).

Overview

Product Name	Anti-PLB1 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Mouse
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q3TTY0

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human PLB1.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PLB1 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.