

Anti-ANKRD6 Antibody

Catalog Number: A08843-1

About ANKRD6

Ankrd6 (Ankyrin repeat domain 6) is a distant homolog of the fly PCP mediator diego and acts as a molecular switch between both Wnt pathways. It promotes signaling of the non canonical Wnt/PCP pathway and simultaneously inhibits gene transcription through the canonical Wnt/beta-catenin pathway.

Overview

Product Name	Anti-ANKRD6 Antibody
Reactive Species	Human
Description	Boster Bio Anti-ANKRD6 Antibody catalog # A08843-1. Tested in WB, IHC, IP, ELISA applications. This antibody reacts with Human.
Application	ELISA, IP, IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9Y2G4

Technical Details

Immunogen	E.coli-derived human ANKRD6 recombinant protein (Position: Q3-L258).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ANKRD6 Antibody

For Research Use Only. Not for use in diagnostic procedures.