

Anti-Peroxin-11B PEX11B Antibody

Catalog Number: A08864-1


Overview

| | |
|----------------------|--|
| Product Name | Anti-Peroxin-11B PEX11B Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-Peroxin-11B PEX11B Antibody catalog # A08864-1. Tested in WB,IHC,ICC,IP applications. This antibody reacts with Human,Mouse,Rat. |
| Application | IP, IHC, ICC, WB |
| Clonality | Polyclonal |
| Formulation | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | O96011 |


Technical Details

| | |
|---------------------|---|
| Immunogen | Synthetic peptide, corresponding to amino acids 132-177 of Human Peroxin-11B. |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | WB: 1:500-1:1000 IHC: 1:50-1:200 ICC: 1:50-1:200 IP: 1:50-1:200 |

Anti-Peroxin-11B PEX11B Antibody (A08864-1) Images



Western blot (WB) analysis of Peroxin-11B polyclonal antibody at 1:500 dilution Lane1:MCF-7 cell lysate Lane2:Raw264.7 cell lysate Lane3:PC12 cell lysate



Immunohistochemistry (IHC) analyzes of Peroxin-11B pAb in paraffin-embedded human liver carcinoma tissue at 1:50 showing Peroxisome membrane staining. Negative control (the right) Using PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG-biotin followed by avidin-peroxidase.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Peroxin-11B PEX11B Antibody

For Research Use Only. Not for use in diagnostic procedures.