

Anti-HNRNPA0 Antibody

Catalog Number: A09015


Overview

Product Name	Anti-HNRNPA0 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-HNRNPA0 Antibody catalog # A09015. Tested in WB,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q13151


Technical Details

Immunogen	Recombinant fusion protein of human HNRNPA0(NP_006796.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IP: 1:50-1:100

Anti-HNRNPA0 Antibody (A09015) Images



Western blot analysis of extracts of various cell lines, using HNRNPA0 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.



Immunoprecipitation analysis of 200ug extracts of 293T cells, using 3 ug HNRNPA0 antibody. Western blot was performed from the immunoprecipitate using HNRNPA0 antibody at a dilution of 1:1000.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HNRNPA0 Antibody

For Research Use Only. Not for use in diagnostic procedures.