

Anti-NREP Antibody

Catalog Number: A09054

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-NREP Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-NREP Antibody catalog # A09054. Tested in WB applications. This antibody reacts with Human,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q16612

Technical Details

Immunogen	A synthetic peptide corresponding to residues in Human NREP.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-NREP Antibody (A09054) Images



Western blot (WB) analysis of NREP polyclonal antibody at 1:500 dilution Lane1:HEK293T whole cell lysate Lane2:H9C2 whole cell lysate

2 Publications Citing This Product

1. PubMed ID: 10.1371/journal.pone.0009995, Investigating the Role of P311 in the Hypertrophic Scar
2. PubMed ID: 2852399, Investigating the Role of P311 in the Hypertrophic Scar

Visit bosterbio.com/anti-nrep-antibody-a09054-boster.html to see all 2 publications.

Anti-NREP Antibody

For Research Use Only. Not for use in diagnostic procedures.