

Anti-Zinc finger homeobox protein 4 ZFHx4 Antibody

Catalog Number: A09122


Overview

Product Name	Anti-Zinc finger homeobox protein 4 ZFHx4 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Zinc finger homeobox protein 4 ZFHx4 Antibody catalog # A09122. Tested in WB, IHC, IF applications. This antibody reacts with Human, Mouse.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q86UP3

Technical Details

Immunogen	Synthesized peptide derived from part region of human protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ZFHx4 antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-1000 IHC 1:200-500 IF 1:50-200

Anti-Zinc finger homeobox protein 4 ZFHx4 Antibody (A09122) Images



Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200 (4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200 (room temperature, 45min).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zinc finger homeobox protein 4 ZFHx4 Antibody

For Research Use Only. Not for use in diagnostic procedures.