

Anti-CCRL2/Cram A Antibody

Catalog Number: A09152


Overview

Product Name	Anti-CCRL2/Cram A Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-CCRL2/Cram A Antibody catalog # A09152. Tested in WB, IF, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O00421

Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human CCRL2. AA range:141-190
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Suggested Dilutions	WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:20000

Anti-CCRL2/Cram A Antibody (A09152) Images



Western Blot (WB) analysis of specific cells using CCRL2 Polyclonal antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CCRL2/Cram A Antibody

For Research Use Only. Not for use in diagnostic procedures.