

Anti-CDK5 Activator-binding C42 CDK5RAP1 Antibody

Catalog Number: A09177-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-CDK5 Activator-binding C42 CDK5RAP1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-CDK5 Activator-binding C42 CDK5RAP1 Antibody catalog # A09177-1. Tested in ELISA, IHC applications. This antibody reacts with Human.
Application	ELISA, IF, IHC
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q96SZ6


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human CDK5RAP1. AA range:251-300
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Suggested Dilutions

IHC 1:100-1:300
ELISA 1:20000
IF 1:50-200

Anti-CDK5 Activator-binding C42 CDK5RAP1 Antibody (A09177-1) Images



Immunohistochemistry validation of CDK5RAP1 using Anti-CDK5 Activator-binding C42 CDK5RAP1 Antibody (A09177-1). Immunohistochemical analysis of paraffin-embedded Human brain. Antibody was diluted at 1:100 (4°C

[Anti-CDK5 Activator-binding C42 CDK5RAP1 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.