

Anti-Helicase POLQ-like HELQ Antibody

Catalog Number: A09185


Overview

Product Name	Anti-Helicase POLQ-like HELQ Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Helicase POLQ-like HELQ Antibody catalog # A09185. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8TDG4

Technical Details

Immunogen	Recombinant fusion protein of human HEL308(NP_598375.2).
Isotype	IgG
Form	Liquid
Concentration	1mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000-1:3000 IHC: 1:50-1:200

Anti-Helicase POLQ-like HELQ Antibody (A09185) Images



Western blot analysis of extracts of various cell lines, using HELQ antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Helicase POLQ-like HELQ Antibody

For Research Use Only. Not for use in diagnostic procedures.