

Anti-BMP-2-inducible protein kinase BMP2K Antibody

Catalog Number: A09187

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-BMP-2-inducible protein kinase BMP2K Antibody
Reactive Species	Human
Description	Boster Bio Anti-BMP-2-inducible protein kinase BMP2K Antibody catalog # A09187. Tested in WB,IP applications. This antibody reacts with Human.
Application	IP, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9NSY1


Technical Details

Immunogen	Synthetic peptide, corresponding to Human BMP2K.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000-1:2000 IP: 1:50-1:200

Anti-BMP-2-inducible protein kinase BMP2K Antibody (A09187) Images



Western blot (WB) analysis of BMP2K polyclonal antibody at 1:1000 dilution
Lane1:A549 whole cell lysate(30ug)
Lane2:MCF-7 whole cell lysate(30ug)



Immunoprecipitation of MCF-7 cell lysates using BMP2K pAb (Sepharose Bead Conjugate) #BD0048 (lane 2 and lane 3) and Nonspecific IgG Control (Sepharose Bead Conjugate) #BD0048 (lane 4). Lane 1 is 30% input. The western blot was probed using BMP2K pAb #BS60064.

Anti-BMP-2-inducible protein kinase BMP2K Antibody

For Research Use Only. Not for use in diagnostic procedures.