

Anti-USP48 Antibody

Catalog Number: A09214


Overview

Product Name	Anti-USP48 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-USP48 Antibody (Catalog# A09214). Tested in IHC-P, IF (paraffin), WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IHC-P, WB, IF (paraffin)
Clonality	Polyclonal
Formulation	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short-term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q86UV5

Technical Details

Immunogen	A peptide derived from human USP48. Immunogen sequence location: 631 - 680
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Suggested Dilutions	Western Blot, 1:500-2000 Immunohistochemistry, 1:100 - 1:300 ELISA, 1:5000 Not yet tested in other applications.

Anti-USP48 Antibody (A09214) Images



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue, using USP48 Antibody. The picture on the right is blocked with the synthesized peptide.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-USP48 Antibody