

Anti-NMUR1/Neuromedin U R1 Antibody

Catalog Number: A09256

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-NMUR1/Neuromedin U R1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-NMUR1/Neuromedin U R1 Antibody catalog # A09256. Tested in WB, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9HB89

Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human NMUR1. AA range:1-50
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml

Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:5000

Anti-NMUR1/Neuromedin U R1 Antibody (A09256) Images



Western Blot (WB) analysis of K562 cells using NMUR1 Polyclonal antibody.

[Anti-NMUR1/Neuromedin U R1 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.