

## Anti-BDH1 Monoclonal Antibody

Catalog Number: A09277


### Overview

Product Name	Anti-BDH1 Monoclonal Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-BDH1 Monoclonal Antibody catalog # A09277. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, IHC, WB
Clonality	Monoclonal 1A5
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q02338


### Technical Details

Immunogen	Purified recombinant fragment of human BDH1 expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Ascitic fluid
Suggested Dilutions	WB 1:500-1:2000 IHC 1:200-1:1000 ELISA 1:10000 IF 1:50-200

## Anti-BDH1 Monoclonal Antibody (A09277) Images



Western Blot analysis using BDH1 Monoclonal Antibody against HepG2 (1) and NIH/3T3 (2) cell lysate.



Immunohistochemistry analysis of paraffin-embedded human liver cancer (left) and colorectal cancer tissues (right) with DAB staining using BDH1 Monoclonal Antibody.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-BDH1 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.