

Anti-Integrin alphaD ITGAD Antibody

Catalog Number: A09283

Overview

Product Name	Anti-Integrin alphaD ITGAD Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Integrin alphaD ITGAD Antibody catalog # A09283. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q13349


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human ITGAD. AA range:901-950
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 IHC: 1:100-1:300 ELISA 1:20000 IF 1:50-200

Anti-Integrin alphaD ITGAD Antibody (A09283) Images



Western blot analysis of lysate from 4T1 cells, using ITGAD Antibody.



Western Blot analysis of 4T1, 293 cells using Integrin alphaD Polyclonal Antibody. Secondary antibody was diluted at 1:20000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Integrin alphaD ITGAD Antibody

For Research Use Only. Not for use in diagnostic procedures.