

## Anti-Phospho-IPP-2 (S120/S121) PPP1R2 Antibody

Catalog Number: A09307S120S121


### Overview

Product Name	Anti-Phospho-IPP-2 (S120/S121) PPP1R2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phospho-IPP-2 (S120/S121) PPP1R2 Antibody catalog # A09307S120S121. Tested in IHC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P41236


### Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human PPP1R2 around the phosphorylation site of Ser120 and Ser121. AA range:86-135
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	IHC 1:100-1:300 ELISA 1:5000 IF 1:50-200

## Anti-Phospho-IPP-2 (S120/S121) PPP1R2 Antibody (A09307S120S121) Images



Immunohistochemistry analysis of paraffin-embedded human heart, using PPP1R2 (Phospho-Ser120+Ser121) Antibody. The picture on the right is blocked with the phospho peptide.



Enzyme-Linked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phosphopeptide (Phospho-left) and Non-Phosphopeptide (Phospho-right), using PPP1R2 (Phospho-Ser120+Ser121) Antibody

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-IPP-2 (S120/S121) PPP1R2 Antibody

For Research Use Only. Not for use in diagnostic procedures.