

## Anti-SAP30L Antibody

Catalog Number: A09316-1

### About SAP30L

SAP30L, also named as NS4ATP2, is one of the key proteins in a multi-subunit protein complex involved in transcriptional regulation via histone deacetylation. SAP30L is an essential component of the Sin3A corepressor complex. The MW migrates to 26-30 kDa with phospho and SUMO modification and may generate complexes.

### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-SAP30L Antibody   |
| Reactive Species     | Human, Mouse, Rat  |
| Description          | Boster Bio Anti-SAP30L Antibody catalog # A09316-1. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.   |
| Application          | ELISA, WB  |
| Clonality            | Polyclonal   |
| Formulation          | 500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol<br>This antibody is supplied in a stabilized formulation.<br>Compatibility with conjugation reactions depends on the chemistry of the conjugation method used.<br>For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | 12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.   |
| Host                 | Rabbit   |
| Uniprot ID           | Q9HAJ7   |

### Technical Details

|                     |  |
|---------------------|--|
| Immunogen           | E.coli-derived human SAP30L recombinant protein (Position: G3-K180). |
| Form                | Liquid   |
| Concentration       | 500 ug/ml  |
| Purification        | Immunogen affinity purified.   |
| Suggested Dilutions | Western blot, 1:500-2000<br>ELISA, 1:100-1000                        |

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-SAP30L Antibody

For Research Use Only. Not for use in diagnostic procedures.