

Anti-GCK/GLK MAP4K2 Antibody

Catalog Number: A09319-1


Overview

Product Name	Anti-GCK/GLK MAP4K2 Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-GCK/GLK MAP4K2 Antibody catalog # A09319-1. Tested in ELISA, WB applications. This antibody reacts with Human, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q12851

Technical Details

Immunogen	Synthesized peptide derived from human GCK/GLK. at AA range: 221-270
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-2000 ELISA 1:10000-20000

Anti-GCK/GLK MAP4K2 Antibody (A09319-1) Images



Western blotting validation for Anti-GCK/GLK MAP4K2
Antibody A09319-1 Western Blot analysis of 1

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GCK/GLK MAP4K2 Antibody

For Research Use Only. Not for use in diagnostic procedures.