

Anti-TMEM119 Antibody

Catalog Number: A09325-1

About TMEM119

Involved in positive regulation of bone mineralization; positive regulation of osteoblast differentiation; and positive regulation of osteoblast proliferation. Located in plasma membrane.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-TMEM119 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-TMEM119 Antibody catalog # A09325-1. Tested in IHC, IF applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, IHC
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q4V9L6

Technical Details

Immunogen	E.coli-derived human TMEM119 recombinant protein (Position: Q120-Q220).
Form	Liquid
Concentration	500 ug/ml

Purification	Immunogen affinity purified.
Suggested Dilutions	Immunohistochemistry, 1:50-400 Immunofluorescence, 1:50-400

[Anti-TMEM119 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.