

Anti-Xaa-Pro aminopeptidase 2 XPNPEP2 Antibody

Catalog Number: A09337

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-Xaa-Pro aminopeptidase 2 XPNPEP2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Xaa-Pro aminopeptidase 2 XPNPEP2 Antibody catalog # A09337. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O43895

Technical Details

Immunogen	Recombinant protein of human XPNPEP2
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000-1:2000

Anti-Xaa-Pro aminopeptidase 2 XPNPEP2 Antibody (A09337) Images



Western blot analysis of extract of various cells, using XPNPEP2 antibody.

Anti-Xaa-Pro aminopeptidase 2 XPNPEP2 Antibody

For Research Use Only. Not for use in diagnostic procedures.