

Anti-HGFA HGFAC Antibody

Catalog Number: A09351

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-HGFA HGFAC Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-HGFA HGFAC Antibody catalog # A09351. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q04756


Technical Details

Immunogen	Synthesized peptide derived from the C-terminal region of human HGFA.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified

Suggested Dilutions

WB 1:500-1:2000
ELISA 1:40000

Anti-HGFA HGFAC Antibody (A09351) Images



Western Blot analysis of extracts from K562, Jurkat cells, using HGFA Polyclonal Antibody. Secondary antibody was diluted at 1:20000

Anti-HGFA HGFAC Antibody

For Research Use Only. Not for use in diagnostic procedures.