

## Anti-Mob1A/B (D127) MOB1B Antibody

Catalog Number: A09405-1


### Overview

Product Name	Anti-Mob1A/B (D127) MOB1B Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Mob1A/B (D127) MOB1B Antibody catalog # A09405-1. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q7L9L4

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 97-157 of Human Mob1A/B.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

## Anti-Mob1A/B (D127) MOB1B Antibody (A09405-1) Images



Western blot (WB) analysis of Mob1A/B (D127) polyclonal antibody at 1:500 dilution Lane1:HeLa whole cell lysate(40ug) Lane2:H1792 whole cell lysate(20ug) Lane3:HCT116 whole cell lysate(40ug) Lane4:CT26 whole cell lysate(40ug) Lane5:PMVEC whole cell lysate(40ug)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Mob1A/B (D127) MOB1B Antibody

For Research Use Only. Not for use in diagnostic procedures.