

Anti-DC-LAMP (E305) LAMP3 Antibody

Catalog Number: A09406-1


Overview

Product Name	Anti-DC-LAMP (E305) LAMP3 Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-DC-LAMP (E305) LAMP3 Antibody catalog # A09406-1. Tested in WB applications. This antibody reacts with Human,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UQV4

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 275-335 of Human DC-LAMP.
Isotype	Rabbit IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000-1:2000

Anti-DC-LAMP (E305) LAMP3 Antibody (A09406-1) Images



Western blot (WB) analysis of DC-LAMP (E305) polyclonal antibody at 1:1000 dilution Lane1:HCT116 whole cell lysate(30ug) Lane2:The Thymus tissue lysate of Rat(40ug) Lane3:A549 whole cell lysate(40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DC-LAMP (E305) LAMP3 Antibody

For Research Use Only. Not for use in diagnostic procedures.