

Anti-IDH3B Antibody

Catalog Number: A09433


Overview

| | |
|----------------------|--|
| Product Name | Anti-IDH3B Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-IDH3B Antibody catalog # A09433. Tested in WB applications. This antibody reacts with Human,Rat,Mouse. |
| Application | WB |
| Clonality | Polyclonal |
| Formulation | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | O43837 |

Technical Details

| | |
|---------------------|---|
| Immunogen | A synthetic peptide corresponding to residues in Human IDH3B |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | WB: 1:500-1:1000 |

Anti-IDH3B Antibody (A09433) Images



Western blot (WB) analysis of IDH3B pAb at 1:5000 dilution
Lane1:H9C2 whole cell lysate(40ug) Lane2:RAW264.7 whole cell lysate(40ug) Lane3:HEK293T whole cell lysate(40ug) Lane4:U-87MG whole cell lysate(40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IDH3B Antibody

For Research Use Only. Not for use in diagnostic procedures.