

## Anti-Na<sup>+</sup> CP type VI alpha SCN7A Antibody

Catalog Number: A09448-1

### Overview

Product Name	Anti-Na <sup>+</sup> CP type VI alpha SCN7A Antibody
Reactive Species	Human, Monkey
Description	Boster Bio Anti-Na <sup>+</sup> CP type VI alpha SCN7A Antibody catalog # A09448-1. Tested in ELISA, IF, WB applications. This antibody reacts with Human, Monkey.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q01118

### Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human SCN7A. AA range:771-820
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-2000 IHC 1:50-300

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Na<sup>+</sup> CP type VI alpha SCN7A Antibody

For Research Use Only. Not for use in diagnostic procedures.