

## Anti-OSBP2 Antibody Picoband®

Catalog Number: A09455-3

### About OSBP2

The protein encoded by this gene contains a pleckstrin homology (PH) domain and an oxysterol-binding region. It binds oxysterols such as 7-ketocholesterol and may inhibit their cytotoxicity. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Overview


|                      |                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-OSBP2 Antibody Picoband®                                                                                                                                                                                                                                                                                                                                                                                                      |
| Reactive Species     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Description          | Boster Bio Anti-OSBP2 Antibody Picoband® catalog # A09455-3. Tested in WB, IHC, ELISA applications. This antibody reacts with Human, Mouse, Rat. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance. |
| Application          | ELISA, IHC, WB                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Clonality            | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Formulation          | Each vial contains 4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na <sub>2</sub> HPO <sub>4</sub> .                                                                                                                                                                                                                                                                                                                                          |
| Storage Instructions | At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.                                                                                                                                                                                                                                  |
| Host                 | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Uniprot ID           | Q969R2                                                                                                                                                                                                                                                                                                                                                                                                                             |

### Technical Details


|               |                                                                           |
|---------------|---------------------------------------------------------------------------|
| Immunogen     | E.coli-derived human OSBP2 recombinant protein (Position: K51-K783).      |
| Form          | Lyophilized                                                               |
| Concentration | Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml. |

|                     |                                                                                                                                              |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------|
| Purification        | Immunogen affinity purified.                                                                                                                 |
| Suggested Dilutions | Western blot, 0.25-0.5 ug/ml, Human, Mouse, Rat<br>Immunohistochemistry(Paraffin-embedded Section), 2-5 ug/ml, Human<br>ELISA, 0.1-0.5 ug/ml |

## Anti-OSBP2 Antibody Picoband® (A09455-3) Images



Western blot analysis of OSBP2 using anti-OSBP2 antibody (A09455-3). Electrophoresis was performed on a 8% SDS-PAGE gel at 80V (Stacking gel) / 120V (Resolving gel) for 2 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: human Hela whole cell lysates, Lane 2: human SiHa whole cell lysates, Lane 3: human K562 whole cell lysates, Lane 4: human U2OS whole cell lysates, Lane 5: rat testis tissue lysates, Lane 6: mouse testis tissue lysates, Lane 7: mouse HEPA1-6 whole cell lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-OSBP2 antigen affinity purified polyclonal antibody (A09455-3) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1% Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an ECL Plus Western Blotting Substrate (Catalog # AR1196-200) with Tanon 5200 system. A specific band was detected for OSBP2 at approximately 101 kDa. The expected band size for OSBP2 is at 101 kDa.



IHC analysis of OSBP2 using anti-OSBP2 antibody (A09455-3). OSBP2 was detected in a paraffin-embedded section of human cerebellum tissue. Heat mediated antigen retrieval was performed in EDTA buffer (pH 8.0, epitope retrieval solution). The tissue section was blocked with 10% goat serum. The tissue section was then incubated with 2 ug/ml rabbit anti-OSBP2 Antibody (A09455-3) overnight at 4°C. Peroxidase Conjugated Goat Anti-rabbit IgG was used as secondary antibody and incubated for 30 minutes at 37°C. The tissue section was developed using HRP Conjugated Rabbit IgG Super Vision Assay Kit (Catalog # SV0002) with DAB as the chromogen.

### Anti-OSBP2 Antibody

For Research Use Only. Not for use in diagnostic procedures.