

## Anti-HDBP1 Antibody (Center)

Catalog Number: A09479

### About SLC2A4RG

THDBP1 is a nuclear transcription factor involved in the activation of the solute carrier family 2 member 4 gene. This protein interacts with another transcription factor, myocyte enhancer factor 2, to activate transcription of this gene.


### Overview

Product Name	Anti-HDBP1 Antibody (Center)
Reactive Species	Human
Description	Boster Bio Anti-HDBP1 Antibody (Center) (Catalog # A09479). Tested in WB, IHC-P application(s). This antibody reacts with Human.
Application	IHC-P, WB
Clonality	Polyclonal
Formulation	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9NR83


### Technical Details

Immunogen	This HDBP1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 213-242 amino acids from the Central region of human HDBP1.
Predicted Reactive Species	Monkey
Isotype	Rabbit IgG
Purification	This antibody is purified through a protein G column, followed by dialysis against PBS.
Suggested Dilutions	WB: 1:1000 IHC-P: 1:10-1:50

## Anti-HDBP1 Antibody (Center) (A09479) Images



HDBP1 Antibody (Center) western blot analysis in SW620 cell line lysates (35ug/lane). This demonstrates the HDBP1 antibody detected the HDBP1 protein (arrow).



Formalin-fixed and paraffin-embedded human lung carcinoma tissue reacted with HDBP1 Antibody (Center), which was peroxidase-conjugated to the secondary antibody, followed by DAB staining. This data demonstrates the use of this antibody for immunohistochemistry; clinical relevance has not been evaluated.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HDBP1 Antibody (Center)

For Research Use Only. Not for use in diagnostic procedures.