

Anti-Guanylate cyclase activator 2B GUCA2B Antibody

Catalog Number: A09489


Overview

Product Name	Anti-Guanylate cyclase activator 2B GUCA2B Antibody
Reactive Species	Mouse
Description	Boster Bio Anti-Guanylate cyclase activator 2B GUCA2B Antibody catalog # A09489. Tested in WB applications. This antibody reacts with Mouse.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q16661

Technical Details

Immunogen	A synthetic peptide of human GUCA2B(NP_009033.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000

Anti-Guanylate cyclase activator 2B GUCA2B Antibody (A09489) Images



Western blot analysis of extracts of mouse small intestine, using GUCA2B antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Guanylate cyclase activator 2B GUCA2B Antibody

For Research Use Only. Not for use in diagnostic procedures.