

Anti-STK36 Antibody

Catalog Number: A09526-2

About STK36

STK36(Serine/threonine-protein kinase 36) is also named as KIAA1278 and belongs to the Ser/Thr protein kinase family. The Serine/threonine protein kinase plays an important role in the sonic hedgehog (Shh) pathway by regulating the activity of GLI transcription factors and controls the activity of the transcriptional regulators GLI1, GLI2 and GLI3 by opposing the effect of SUFU and promoting their nuclear localization. It has 2 isoforms produced by alternative splicing.

Overview

Product Name	Anti-STK36 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-STK36 Antibody catalog # A09526-2. Tested in WB, IHC, IP, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IP, IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9NRP7

Technical Details

Immunogen	E.coli-derived human STK36 recombinant protein (Position: H5-I301).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 ImmunoPrecipitation, 1:250-300 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-STK36 Antibody

For Research Use Only. Not for use in diagnostic procedures.