

Anti-ADGRA2 Antibody Picoband®

Catalog Number: A09540-2

About ADGRA2

Predicted to enable G protein-coupled receptor activity. Involved in canonical Wnt signaling pathway. Located in plasma membrane. Part of Wnt signalosome.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-ADGRA2 Antibody Picoband®
Reactive Species	Rat
Description	Boster Bio Anti-ADGRA2 Antibody Picoband® catalog # A09540-2. Tested in WB applications. This antibody reacts with Rat. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	WB
Clonality	Polyclonal
Formulation	Each vial contains 4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q96PE1

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human ADGRA2.(DMMDVVYVAQMIQKFL).
Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.

Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.25-0.5 ug/ml, Rat

[Anti-ADGRA2 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.