

## Anti-IQGAP3 Antibody

Catalog Number: A09547

### Overview


Product Name	Anti-IQGAP3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-IQGAP3 Antibody catalog # A09547. Tested in WB, ELISA, IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q86VI3

### Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human IQGAP3. AA range:107-156
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Suggested Dilutions	WB 1:500-2000 IHC 1:50-300 ELISA 1:2000-20000

## Anti-IQGAP3 Antibody (A09547) Images

---



Western blot analysis of various cells using IQGAP3 Polyclonal Antibody diluted at 1:1000

### Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-IQGAP3 Antibody

For Research Use Only. Not for use in diagnostic procedures.