

## Anti-FPR3 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A09576-1-Fluoro594

### About FPR3

Predicted to enable N-formyl peptide receptor activity and complement receptor activity. Predicted to be involved in several processes, including complement receptor mediated signaling pathway; phospholipase C-activating G protein-coupled receptor signaling pathway; and positive regulation of cytosolic calcium ion concentration. Predicted to act upstream of or within G protein-coupled receptor signaling pathway. Predicted to be located in membrane. Predicted to be active in plasma membrane.

### Overview

Product Name	Anti-FPR3 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P25089

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human FPR3.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-FPR3 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.